






# Biocatalyst Design








## Overlap extension primer



**Figure 2.** Overlap extension PCR method: → represents a primer, and × represents a mutagenic codon.


## Whole plasmid single round PCR




Two oligonucleotide primers containing the desired mutation(s), X, are used with a DNA polymerase to extend the original plasmid in single PCR reaction, resulting in a nicked circular strand.

DpnI methylase digests the methylated, nonmutated parental DNA template.

The circular, nicked dsDNA is transformed into competent cells. The competent cells repair the nicks in the mutated plasmids.



## Whole Gene: mutation





-> PCR with low fidelity !!!

Achieved by:


- Increased  $Mg^{2+}$  concentration
- Addition of  $Mn^{2+}$
- Not equal concentration of the four dNTPs
- Use of dITP
- Increasing amount of Taq polymerase (Polymerase with NO proof reading function)

# Single AA/site mutation: degenerative primers





## DNA shuffling




DNase I treatment (Fragmentation,  
10-50 bp,  $Mn^{2+}$ )

Reassembly (PCR without primers,  
Extension and Recombination)

PCR amplification

## Family shuffling


Native genes



DNA fragmentation  
and PCR



Hybrid genes




Genes coming from the same  
gene family -> highly  
homologous

-> Family shuffling

# Directed Evolution

Difference between non-recombinative and recombinative methods



Non-recombinative methods


recombinative methods ->  
hybrids (chimeric proteins)

# Screening:

## Basis for all screening & selection methods

Expression Libraries

->link gene with encoded product which is responsible for enzymatic activity



# Low-medium throughput screens

- > Detection of enzymatic activity of colonies on agar plates or "crude cell lysates" -> production of fluorophor or chromophor or halos
- > Screen up to  $10^4$  colonies
- > effective for isolation of enzymes with improved properties
- > not so effective for isolation of variants with dramatic changes of phenotype






Lipase: variants on Olive oil plates  
With pH indicator (brilliant green)

## Phage display screening


### Filamentous phage display

Bacteriophage surface display is a system for the expression of (poly-)peptides (phenotype) that are physically linked to their coding DNA (genotype). This is achieved by cloning a foreign gene in fusion to a phage coat protein, usually pIII or pVIII (below).




# HTS of enzymes with Cell Display

- **Yeast display** of antibody scFv fragments
  - > Expression of scFv on yeast cells is monitored using either the HA or c-myc epitope tags.
  - > Binding of the target antigen, HIV-1 gp120, is visualized using a biotinylated mAb to a noncompetitive epitope on gp120 and fluorescent streptavidin. Gp120-binding scFvs are selected by fluorescence activated cell sorting (FACS) of the yeast cells.










# In vitro compartmentalization (IVC) or Cell free systems

- Water-in-oil emulsion -> make microscopic compartments (droplet volume ~5 fL) -> no diffusion between compartments.
- Each compartment contains in general one single gene and acts as an artificial cell (invitro transcription and translation)
- -> Gene is linked to substrate
- -> if protein is active substrate will be converted into product -> signal -> fishing out
- -> direct access to gene




# In vitro compartmentalization (IVC)


-> w/o/w -> water in oil in water emulsion (Vesicle)

Avantage:

- No direct link between product and gene necessary -> keep the compartment
- Can be directly analysed with FACS



# FACS (Fluorescens-activated cell sorter)



Capacity >  $10^7$  per hour

- The cell suspension is entrained in the center of a narrow, rapidly flowing stream of [liquid](#).
- A [vibrating](#) mechanism causes the stream of cells to break into individual droplets.
- Just before the stream breaks into droplets, the flow passes through a fluorescence measuring station where the fluorescent character of interest of each cell is measured.
- An electrical charging ring is placed just at the point where the stream breaks into droplets. A [charge](#) is placed on the ring based on the immediately-prior fluorescence intensity measurement, and the opposite charge is trapped on the droplet as it breaks from the stream. The charged droplets then fall through an [electrostatic deflection](#) system that diverts droplets into containers based upon their charge.



# Lab on a Chip

PAPER


[View Article Online](#)

[View Journal](#) | [View Issue](#)


## High-throughput screening for industrial enzyme production hosts by droplet microfluidics†

Cite this: *Lab Chip*, 2014, 14, 806


Staffan L. Sjoström,<sup>a</sup> Yunpeng Bai,<sup>a</sup> Mingtao Huang,<sup>b</sup> Zihé Liu,<sup>b</sup> Jens Nielsen,<sup>abc</sup>  
Haakan N. Joensson<sup>a</sup> and Helene Andersson Svahn<sup>\*a</sup>




b)







c)



d)





**a)****b)****c)**



