



# ASSAY OF B- GLUCOSIDASE

---

LAXMAN MANJHI 2019BBI0034

RATNESH KUMAR SHARMA 2019BBI0047



# BACKGROUND

---

- $\beta$ -Glucosidase is an enzyme that catalyzes the hydrolysis of  $\beta$ -D-glucosides into glucose and another molecule. An assay of  $\beta$ -glucosidase measures the amount of enzyme activity in a sample. This is typically done by measuring the amount of glucose produced from a substrate over time and converting that measurement to units of activity. The assay can be used to determine the activity of  $\beta$ -glucosidase in various samples, including bacteria, yeast, plant tissues, and more. The specific conditions of the assay, such as temperature, pH, and substrate concentration, may be optimized for each sample to ensure accurate results.

# METHODS

---


- I. Standard Curve a. Dilute PNP solution from 10  $\mu\text{M}$  to 100  $\mu\text{M}$  with citrate buffer to a final volume of 1 mL. You can use dilutions of 20, 40, 60, 80, or 100  $\mu\text{M}$ . (Note: Dilute with citrate buffer, not water) b. Add 0.5 mL of 1N  $\text{Na}_2\text{CO}_3$  to each tube and measure the OD at 405 nm against a reagent blank (Blank = 1 mL citrate buffer + 0.5 mL 1N  $\text{Na}_2\text{CO}_3$ ).
  
- II. Enzyme Assay a. Add 0.9 mL of pre-equilibrated PNPG solution to 100  $\mu\text{L}$  of enzyme. b. Incubate the solution in a water bath for 10 minutes. c. Stop the reaction by adding 0.5 mL of 1N  $\text{Na}_2\text{CO}_3$  to each tube. d. Measure the absorbance at 405 nm. (Blank = 0.1 mL citrate buffer + 0.9 mL PNPG + 0.5 mL of 1N  $\text{Na}_2\text{CO}_3$ )

# OBSERVATION

---

conc	abs	
	0	0
	20	0.27
	40	0.535
	60	0.725
	80	1.047
	100	1.279
Unknown		
50c		1.154
RM		0.443

# GRAPH



# CALCULATION

---

Equation is  $y = 0.0128x$

Where,

$Y = \text{Abs}$

$X = \text{Conc.}$

Conc. At 50c = 90

Conc. At RT = 34.60

---

THANK YOU