



est strips w/ Long shelf-life of upto
ifetime Warranty



Enzyme Biosensors

Ravikrishnan Elangovan,
Department of Biochemical Engg and Biotechnology
Indian Institute of Technology - Delhi




Brief history of Biosensors



- 1916 First report on immobilization of proteins : adsorption of invertase on activated charcoal
- 1922 First glass pH electrode
- 1956 Clark published his definitive paper on the oxygen electrode.
- 1962 First description of a biosensor: an amperometric enzyme electrode for glucose (Clark)
- 1969 Guilbault and Montalvo – First potentiometric biosensor: urease immobilized on an ammonia electrode to detect urea
- 1970 Bergveld – ion selective Field Effect Transistor (ISFET)
- 1975 Lubbers and Opitz described a fibre-optic sensor with immobilised indicator to measure carbon dioxide or oxygen.

Structure of a Biosensor



A biosensor is an analytical device which is used to determine the presence and concentration of a specific substance in a biological analyte

<http://dx.doi.org/10.15226/2471-3627/4/1/00116>

What are different types of Biomarker ?






Samples :

1. Blood
2. Urine
3. Fecal
4. Nasal swab
5. Body fluids
6. Tissue samples
7. Environmental sample
8. Breath
9. Food


What could be biomarkers??

1. Protein
2. DNA/RNA
3. Lipids
4. Metabolites
5. Ions/metals


Recognizing element



Bioreceptor: Antibody



Monoclonal Antibody Binding



Polyclonal Antibody Binding




1 Inject antigen into rabbit.


2 Antigen activates B cells.

3 Plasma B cells produce polyclonal antibodies.

4 Obtain antiserum from rabbit containing polyclonal antibodies.





Bioreceptor: Monoclonal Antibody



<https://microbenotes.com/monoclonal-antibodies-types-uses-and-limitations/>

Bioreceptor: Aptamers

A



<https://doi.org/10.1523/JNEUROSCI.1969-16.2017>

<https://doi.org/10.1007/s00216-007-1346-4>

Course evaluation plan

So far....

Minor 1 - 20 marks
Quiz - 10 marks
Lab components - 10 marks

For graduating students

Do it from home Minor 2 – 15 marks

Oral viva - 10 marks

Lab notebook - 10 marks

Business plan/Term paper – 25 marks

Lecture plan:


Biosensor complete – 13th June
Bio fuels complete - 20th June (2 lectures)
Techno-commercial – 27th June (2 lectures)

All lecture components will be completed before 30th June!


Request all the graduating students to complete submission by 5th July

Request all the others students to complete submission by 15th July

Transducer



Optical transducer: Absorption





$$I_1/I_0 = e^{-\alpha lc}$$

l is the pass length


C is the concentration of absorbing material

α is the absorption coefficient




Color reaction



Optical transducer: Fluorescence



An optical system for
Fluorescence measurement



Optical transducer: Fluorescence



© 2004 Thomson - Brooks/Cole

Optical transducer: Fluorescence probes




Fig. 3. (a) Ethidium bromide; (b) the process of intercalation, illustrating the lengthening and untwisting of the DNA helix.




Figure 10.2.15 Principle of coupled enzymatic assays using our Amplex® Red reagent. Oxidation of glucose by glucose oxidase results in generation of H_2O_2 , which is coupled to conversion of the Amplex® Red reagent to fluorescent resorufin by HRP. The detection scheme shown here is used in our Amplex® Red Glucose/Glucose Oxidase Assay Kit (A22189).

Optical transducer: Fluorescence probes






Figure 10.1.16 Principle of the fluorogenic response to protease cleavage exhibited by HIV protease substrate 1 (H2930). Quenching of the EDANS fluorophore (F) by distance-dependent resonance energy transfer to the dabcyI quencher (Q) is eliminated upon cleavage of the intervening peptide linker.



FRET pairs: Peptide



FRET pairs: Taqman probe


Reference: Thermofisher.com

ELISA



<https://www.abcam.com/kits/elisa-principle>

Types of ELISA



COVID Virus Dx test
Simultaneous absorbion



COVID Antibody test

COVID Virus Dx test

COVID Virus Dx test,
Sites have competitor molecules



Basic principles of ELISA

Sandwich ELISA



- the antibodies fixed to a solid surface, such as the inner surface of a test tube;
- a preparation of the same antibodies coupled to an enzyme.

Amplification effect of Enzyme in biosensing



High sensitivity: $K_{cat} > 100 \text{ s}^{-1}$: results are amplified!
High specificity: Enzyme have high selectivity





Point-of-care diagnostics for noncommunicable diseases using synthetic urinary biomarkers and paper microfluidics

Andrew D. Warren^{a,b,1}, Gabriel A. Kwong^{a,b,1}, David K. Wood^{c,1}, Kevin Y. Lin^{b,d}, and Sangeeta N. Bhatia^{a,b,e,f,g,h,2}


^aHarvard–Massachusetts Institute of Technology Division of Health Sciences and Technology, Institute for Medical Engineering and Science, Massachusetts Institute of Technology, Cambridge, MA 02139; ^bDavid H. Koch Institute for Integrative Cancer Research, Massachusetts Institute of Technology, Cambridge, MA 02139; ^cDepartment of Biomedical Engineering, University of Minnesota, Minneapolis, MN 55455; ^dChemical Engineering, Massachusetts Institute of Technology, Cambridge, MA 02139; ^eElectrical Engineering and Computer Science, Massachusetts Institute of Technology, Cambridge, MA 02139; ^fDepartment of Medicine, Brigham and Women's Hospital and Harvard Medical School, Boston, MA 02115; ^gBroad Institute of Massachusetts Institute of Technology and Harvard, Cambridge, MA 02139; and ^hHoward Hughes Medical Institute, Cambridge, MA 02139

Edited by Stephen R. Quake, Stanford University, Stanford, CA, and approved January 21, 2014 (received for review August 1, 2013)

Case study 1: Enzyme sensor for cancer detection



Case study 1: Enzyme sensor for cancer detection



ARTICLES


PUBLISHED ONLINE: 2 SEPTEMBER 2012 | DOI: 10.1038/NCHEM.1435

nature
chemistry


Rapid point-of-care detection of the tuberculosis pathogen using a BlaC-specific fluorogenic

Hexin Xie^{1†}, Joseph Mire^{2†}, Ying Kong^{3†}, MiHee Chang³, Hany A. Hassounah³, Chris N. T. James C. Sacchettini², Jeffrey D. Cirillo³ and Jianghong Rao^{1*}

These probes showed an enhancement by 100–200 times in fluorescence emission on BlaC activation and a greater than 1,000-fold selectivity for BlaC over TEM-1 β -lactamase



CDC-1: $R^1 = \text{Bn}$, $R^2 = \text{H}$, $n = 0$
CDC-2: $R^1 = \text{Bn}$, $R^2 = \text{H}$, $n = 1$




CDC-3: $R^1 = \text{EtN}$, $R^2 = \text{H}$, $n = 0$

CDC-4: $R^1 = \text{Me}$, $R^2 = \text{H}$, $n = 0$
CDC-5: $R^1 = \text{Me}$, $R^2 = \text{OMe}$, $n = 0$
CDC-OMe: $R^1 = \text{Bn}$, $R^2 = \text{OMe}$, $n = 0$
CDC-OEt: $R^1 = \text{Bn}$, $R^2 = \text{OEt}$, $n = 0$

Figure 1 | General structures of blue fluorescent probes and their hydrolysis by β -lactamase (Bla), which triggers the release of umbelliferone and turns on fluorescence.

Case study 2: Enzyme sensor for TB detection





Thank you